

Certificate of Compliance


SOTAX

Solutions for Pharmaceutical Testing

SOTAX AG, Nordring 1, Aesch 4147, Switzerland

Part: **1 Liter Vessel**
Number: 13411-01
Serial Number: 74641 to 74650

For SOTAX AT Dissolution Testing Systems



Dimensions:	Inside Ø	a	100 mm	± 0.3 mm
	Height	b	186.9 mm	± 0.7 mm

With hemispherical bottom

Material: Borosilicate glass

SOTAX Quality Management
Bruno Balli
2018-05-24


Signature